# KAnalyze Manual v0.9.7

# Contents

1	Intr	oducti	on	3							
	1.1	About	K-mers	3							
	1.2	About	KAnalyze	3							
<b>2</b>	Command Line Usage										
	2.1	Runnii	ng Modules	4							
	2.2	Java C	ommand	4							
	2.3	Count	Module	4							
		2.3.1	Usage	4							
		2.3.2	Example Usage	6							
		2.3.3	Properties	6							
		2.3.4	Segment Size Format	6							
		2.3.5	Count Performance Tuning	7							
	2.4	Stream	ı Module	8							
		2.4.1	Usage	8							
		2.4.2	Example Usage	8							
		2.4.3	Properties	8							
	2.5		Interpolation of the control of the	9							
	2.0	2.5.1	Usage	9							
		2.5.2	Example Usage	9							
			2	Ü							
3	$\mathbf{AP}$			10							
	3.1		oc Documentation	10							
	3.2		t Organization	10							
	3.3	Compo	onents	11							
		3.3.1	Count Merge Component	11							
		3.3.2	Count Split Component	11							
		3.3.3	Count File Writer	12							
		3.3.4	File Writer	12							
		3.3.5	K-mer	12							
		3.3.6	Reverse Complement	13							
		3.3.7	Reader	13							
	3.4	Utility	Classes	14							
		3.4.1	BoundedQueue	14							
		3.4.2	KmerUtil	14							
		3.4.3	LongBoundedQueue	14							
		3.4.4	LongCounter	14							
		3.4.5	LongHashSet	14							
		3.4.6	LongIterator	15							
		3.4.7	StringUtil	15							
	ъ.	7.		1.0							
4			KAnalyze	16							
	4.1		es	16							
	4.2	Sequen	ace Readers	17							

<b>5</b>	$\mathbf{Sup}$	plementary Information 1	18
	5.1	Input File Formats	18
		5.1.1 FASTA	18
		5.1.2 FASTAGZ	18
		5.1.3 FASTQ	18
		5.1.4 FASTQGZ	18
		5.1.5 RAW	18
	5.2	Kmer Output Formats	19
		5.2.1 Sequence Format	19
			19
		5.2.3 Hexadecimal Format	19
		5.2.4 FASTA Format	19
	5.3	Command Line Return Codes	19
	5.4	Building From Source	20
			20
	5.5	Running Unit Tests	21
6	Lice	ense 2	22
	6.1	Documetation	22
	6.2	KAnalyze Software	22

# 1 Introduction

#### 1.1 About K-mers

K-mers are short nucleotide sequences of a fixed length, k. For example, short sequences of length 8 are 8-mers.

Converting sequences to k-mers is an iterative process. To get the first k-mer, start at the first base and read the first k characters. The second k-mer starts at the second base, and so on until the end of the last k-mer reaches the end of the sequence.

Analyzing sequence data sets as k-mer sets has many applications in bioinformatics. Some sequence assembly and gene annotation algorithms rely on k-mers.

# 1.2 About KAnalyze

KAnalyze converts nucleotide (DNA and RNA) sequences to k-mers. KAnalyze evaluates 5 different bases, A, C, G, T, and U. All other letters are skipped and omitted from any k-mer. It can be run as a command line utility, or integrated into pipelines through its command line interface or programming interface (API).

This program was built for speed, accuracy, and longevity. It can compete with the fastest k-mer programs available, it has been extensively tested, and the source code is built to survive updates.

KAnalyze is written in Java, and it requires Java 7 to run. Any machine with a Java 7 JVM (Java Virtual Machine) should be able to run KAnalyze. We have tested it on Linux, Windows, and MacOS.

KAnalyze can currently generate k-mers of size 1 to 31. To optimize for speed, k-mers are processed numerically, which creates this 31-mer size restriction.

# 2 Command Line Usage

KAnalyze has several modules, and each module executes a different task. Each module is a complete pipeline that takes raw input and writes the final results. For example, the *count* module reads sequence files and writes a sorted tab-delimited file of all k-mers and their counts. This section covers how to run each module.

Java 1.7 (aka Java 7) must be installed in order to run any commands in this section. Run java -version to see your version of Java. If you do not see "java 1.7.0" or a later version, Java must be updated before continuing.

# 2.1 Running Modules

KAnalayze is packaged with a program to run each module. The Windows package comes with exe files, and the Linux package comes with Bash scripts. These programs setup the default maximum memory for the program and runs the module.

For example, the count module can be run by typing 'count' from the command line. On Linux, you may need to begin the command with './'1 ('./count' in this case).

# 2.2 Java Command

In most cases, using the module programs described in Section 2.1 is adequate. Sometimes, it's necessary to run the jar file directly.

The current built-in modules allocate 2GB memory. Some advanced parameters, such as multi-threading components or increasing queue sizes, requires more than 2GB of memory. In these cases, the module programs will not work and the jar file must be run directly.

General form:

```
java -jar -Xmx2G kanalyze.jar <module> [<module arguments>]
```

This allocates 2GB of RAM (-Xmx2G), which is required by the built-in modules. More memory can be allowed by increasing 2G.

After the <module> argument, the rest of the command line is the same as above.

# 2.3 Count Module

The count module counts k-mers from all input files. This module is suitable for large sequence data sets. This module uses a fixed amount of memory to accumulate k-mers. When that memory is full, k-mers are sorted and k-mer counts are written to a file. The memory is refilled with the next set of k-mers, and the process repeats. After all k-mers have been written to files, each file is opened and the k-mer counts are merged and written to the output file. This module is designed to use about 2GB of memory.

#### 2.3.1 Usage

```
count [-d <k-mer threads>] [-f <format>] [-g <segsize>] [-k <ksize>]
        [-l <split threads>] [-m <outfmt>] [-o <out>] [-p <key=value>]
        [-r <reverse>] [-r] [-R] [-s] [-S] [-t <temploc>] [-v] [-V] [-x] [-X]
count -h
-d --kmerthreads [default = 1]
```

Set the number of threads to use for the pipeline k-mer step.

 $<sup>^{1}</sup>$ It depends on how your PATH environment variable is setup. If the current directory (.) is in your path, then './' is not necessary.

#### -f --format [default = raw]

Set the input sequence format type. This option determines how the format files are read. Valid examples include "raw", "fasta", "fastagz", and "fastqgz". This option may be set multiple times when reading multiple files with different formats.

#### -g --segsize [default = 65536000]

Size of k-mer segments stored in memory. This many k-mers are accumulated into an array in memory. When the array is full, it is sorted, written to disk, and cleared for the next set of k-mers. This is a performance tuning parameter that does not need to be set for the majority of data sets. If a value is chosen that is higher than the default, then more than 2GB of memory will be required. See the KAnalyze manual for more information on performance tuning. This number may be in decimal, hexadecimal, or octal. An optional multiplier (k, m, or g) may included after the number. For the maximum value, use keyword "max". See the KAnalyze manual usage documentation for more information on acceptable formats.

#### -h --help

Print help and exit.

## -k --ksize [default = 31, max = 31]

K-mer size. This determines the size of k-mers extracted from input sequences.

## -l --splitthreads [default = 1]

Set the number of threads to use for the pipeline split step.

#### -m --outfmt [default = seq]

Sets the format of the output file. Several output formats are built into KAnalyze, including seq, hex, dec, and fasta.

# -o --out [default = count.kc]

Output file name for the file writer.

# -p --property

Set a key/value pair for the modules.

#### -r --reverse

Reverse complement k-mers as they are generated. This option has an optional mode argument. Valid modes are "duplicate", "single", and "canonical". In the default mode, duplicate, the original k-mers and their reverse complements are counted. Mode single counts only reverse complements. Mode canonical counts either the original k-mer or its complement, whichever comes first in ASCII sort order. When specifying a mode, do not include a space after '-r' (e.g. -rcanonical).

#### -R --noreverse [default]

Disable reverse complementing k-mers as they are counted.

#### -s --splitonly

Read files, split k-mers into sorted segments, but do not merge. This option is useful for creating split segments on a distributed system, but merging must take place on the same machine. For more information on splitting and merging in separate steps, see the count module usage information in the KAnalyze manual.

# -S --nosplitonly [default]

Perform splitting and merging in the same operation. This option negates  $\neg s$ 

## -t --temploc [default = Output file directory]

The location where segments are offloaded. This argument must be a directory or the location for a new directory. Parent directories will be created as needed.

#### -v --verbose

Output extra information after job is complete.

#### -V --noverbose [default]

Disable extra information output.

## -x --autodelete [default]

Automatically delete segment (temporary) files generated by the split component after merging. In -s (--splitonly) mode, segment files are not deleted. Segment files input on the command line are also not automatically deleted.

#### -X --noautodelete

Do not automatically delete segment files on the command line.

## 2.3.2 Example Usage

java -jar kanalyze.jar count -k 8 -o counts.kc -f fasta test.fa

Reads sequences from test.fa, counts all 8-mers, and writes k-mer counts to counts.kc.

java -jar kanalyze.jar count -k 8 -o counts.kc -m hex -f fasta test.fa

Reads sequences from test.fa, counts all 8-mers, and writes k-mer counts to counts.kc as hexadecimal strings.

## 2.3.3 Properties

This module recognizes the following properties (-p)

Property	Default	Description
comp.count.autodelete	true	Automatically delete segment files generated by the split
		component. This does not affect segment files input on the
		command line. This parameter is automatically set by -x
		and -X
comp.count.segsize	98304000	Count split component writes segments of this many k-
		mers. Increasing this parameter requires additional mem-
		ory. This parameter is automatically set by -g
comp.kmer.batchsize	100000	The size of k-mer batches passed to the next component
		after k-merizing reads.
kanalyze.outfmt	seq	K-mer output format. Valid values are "seq", "dec", "hex",
		and "fasta". Setting this property has the same effect as
		the -m (-outfmt) option.

# 2.3.4 Segment Size Format

The segment size argument (-s or --segsize) is an integer number represented in decimal, hexadecimal, or octal and optionally followed by a multiplier. Acceptable numbers begin with an optional sign (+ or -) followed by a decimal, hexadecimal, or octal representation of a number. Hexadecimal numbers must begin with "0x", "0X", or "#". Octal numbers must begin with a zero.

Valid multipliers are:

• k or kb: Binary kilo (1,024)

• kd: Decimal kilo (1,000)

• m or mb: Binary mega (1,024<sup>2</sup> or 1,048,576)

• md: Decimal mega (1000000)

• g or gb: Binary giga (1024<sup>3</sup> or 1,073,741,824)

• gd: Decimal giga (1,000,000,000)

The keyword "max" may be used as the whole argument. The maximum value is 2,147,483,647<sup>2</sup>.

No part of this number string is case sensitive (integer or multiplier).

#### 2.3.5 Count Performance Tuning

The count module has three performance tuning parameters. They are -d (--kmerthreads), -l (--splitthreads), and -s (--segsize). Use of these parameters will not affect the final results, but it will affect resource usage (memory, CPU, and I/O) the amount of time required complete the task. The default values for these parameters work on most systems and is suitable for most data sets. When performance tuning KAnalyze, it is recommended to do a default run on an anticipated data set first, then try the tuning parameters and compare the time required to complete the task.

KAnalyze operates in several steps. They are:

- 1. Extract sequence reads from source files.
- 2. Convert sequence reads to k-mers.
- 3. Split k-mers into sorted segments.
- 4. Merge all sorted segments.
- 5. Write sorted k-mer counts to disk.

The -d (--kmerthreads) parameter controls how many threads are used to convert sequence reads to k-mers. The conversion algorithm is very fast, and this is often not main performance bottle-neck. Increasing the number of threads may reduce the amount of time the split component waits for k-mers. If the number of split threads is increased (see below), then increasing this parameter may improve performance.

The -l (--splitthreads) parameter controls how many threads are used to accumulate and write sorted k-mer segments to disk. The split component takes k-mers from the k-mer component, stores them in memory, and when the amount of memory allocated to k-mers is full, it writes the sorted segment to disk. The time required for this component to accumulate k-mers, sort them, and perform the required disk I/O can stop the KAnalyze process while it completes. Having two or more split components may increase performance so one can be accumulating an sorting while another is writing. Too many threads may have a negative performance impact by writing to multiple files at once or by excessive context switching (as active threads change).

The -s (--segsize) parameter controls how many k-mers are accumulated in memory by the split component before they are sorted and offloaded to disk. The parameter argument is the size of a k-mer array. When that array is full, it is sorted and written. The actual memory taken by this array is this argument multiplied by 8 (k-mers are long integers and require 8 bytes each). The maximum value ("-s max") will allocate 16GB of memory to the array. On very large data sets, the default segment size may result in an excessive number of files, which slows down the merging process. In that case, increasing this value may improve performance. If the memory usage is not adequately increased with this parameter, excessive GC (Garbage Collection) activity will degrade performance or terminate the operation.

Increasing any of these parameters from their default values requires additional memory above the 2GB default limit. See Section 2.2 for more information on how to run KAnalyze with additional memory

<sup>&</sup>lt;sup>2</sup>Java's array size limit.

allocated for the JVM (Java Virtual Machine). When the memory limitation is reached, "Out of Memory" or "GC overhead limit exceeded" errors occur. Note that "GC" in this error message refers to "Garbage Collection" and not the nucleotide bases "G" and "C"<sup>3</sup>.

#### Stream Module

The stream module converts sequences to k-mers and dumps them to a file. This is most likely used to generate k-mers for another program to analyze.

#### 2.4.1 Usage

```
stream [-f <format>] [-k <kmer size>] [-m <outfmt>] [-o <output>]
       [-p <key=value>] -v -V
stream -h
-f --format [default = raw]
        Set the input sequence format type. This option determines how the
        format files are read. Valid examples include "raw", "fasta",
        "fastq", "fastagz", and "fastqgz". This option may be set multiple
        times when reading multiple files with different formats.
-h --help
        Print help and exit.
-k --ksize [default = 31, max = 31]
        K-mer size. This determines the size of k-mers extracted from
        input sequences.
-m --outfmt [default = seq]
        Sets the format of the output file. Several output formats are built
        into KAnalyze, including seq, hex, dec, and fasta.
-o --out
```

Output file name for the file writer.

-p --property

Set a property in the form "key=value".

-v --verbose

Output extra information after job is complete.

-V --noverbose [default]

Disable extra information output.

# 2.4.2 Example Usage

java -jar kanalyze.jar stream -k 8 -o stream.km -f fasta test.fa Converts test.fa to 8-mers and writes each 8-mer to stream.km

#### 2.4.3 Properties

This module recognizes the following properties (-p)

 $<sup>^3</sup>$ Unfortunately, this is a JVM-generated error message we have no control over.

Property	Default	Description
comp.count.deletemp	true	Count merge component deletes temporary files after com-
		pleting.
comp.count.segsize	98304000	Count split component writes segments of this many k-
		mers. Increasing this parameter requires additional mem-
		ory.
kanalyze.outfmt	seq	K-mer output format. Valid values are "seq", "dec", "hex",
		and "fasta".
comp.kmer.batchsize	100000	The size of k-mer batches passed to the next component
		after k-merizing reads.

# 2.5 Cite Module

Outputs information on citing KAnalyze.

# 2.5.1 Usage

cite [ref|bibtex]
cite -h

-h --help

Print help and exit.

Outputs information on citing this software in papers and publications.

Valid citation types are: ref, bibtex

# 2.5.2 Example Usage

java -jar kanalyze.jar cite

# 3 API

KAnalyze is not only a command line program, it is also an API (Application Programming Interface) that can be driven from other programs. The CLI (Command Line Interface) itself is a simple application that calls the API to carry out tasks

KAnalyze is a Java program, so the API is implemented as a set of Java classes with well documented interfaces. Any Java program can directly call the API classes. Other programs can use technologies such as JNI (Java Native Interfaces) to drive the API from native code, such as C or C++.

Lastly, the command line itself can be used to drive the program from a script or batch file. Because KAnalyze always terminates with a well-defined return code, scripts can easily execute the CLI and check for errors. For a list of the return codes, see section 5.3.

The remainder of this section outlines the structure and capabilities of the API. It is assumed that a reader of this section has at least a fundamental understanding of Java.

## 3.1 Javadoc Documentation

The KAnalyze API is fully documented with Javadoc comments. Every class and class member (method and field), regardless of scope, has a Javadoc comment. All method arguments, return values, and exceptions are contained in the method's Javadoc comment. For any method arguments that are objects, the comment must disclose what happens when null values are received for that field (either by the argument's comment, or the comment on the exception if one is thrown). All exceptions, whether or not they are runtime exceptions, must be documented along with the conditions that cause them to be thrown. Assertion errors are not documented. Any deviations from these rules should be reported as a program bug.

Javadoc comments can be converted into a set of HTML pages that document the code. In fact, the Java API itself uses Javadoc, so most Java programmers are familiar with the format of the HTML pages. The KAnalyze build system will create the Javadoc pages for two different levels: api and full. The API level contains only public and protected members, while the full documentation contains everything. For a programmer interested in using the API or creating a custom component, the API level documentation is probably sufficient. The full documentation is intended for maintenance programmers.

For more information on building Javadoc HTML pages from source, see section 5.4.1.

## 3.2 Project Organization

The KAnalyze project is organized into several packages, which are all sub-packages of edu.gatech.kanalyze. The major packages are discussed in this section. Most of these major packages also have sub-packages.

The main packages are:

edu.gatech.kanalyze.comp Components that are used to implement in-memory pipelines. Most pipeline components take input from one component, process it, and output to another component. Linking them together forms a pipeline in memory where each pipeline has its own thread. The bulk of the work performed by KAnalyze is done by these components.

edu.gatech.kanalyze.io IO classes used by multiple parts of the KAnalyze system.

edu.gatech.kanalyze.module Modules form the major commands that are run on the command line. Modules load and link components in a pre-defined way. They can be thought of as complete pipelines. The command line interface runs these modules.

edu.gatech.kanalyze.test JUnit test packages. To run tests, see section 5.5.

edu.gatech.kanalyze.util General purpose utility classes run by many parts of the KAnalyze system.

# 3.3 Components

Components are the pipeline elements that are connected together. They are designed to read and write from synchronized queues. Each component instance should run in its own thread. Because of their parallel nature, downstream components can process data from upstream components as it becomes available. The pipeline architecture is the key to the flexibility and speed of KAnalyze.

Reading and writing individual elements from synchronized queues as they are passed from one element to another is inefficient. The function calls are synchronized, so repetitive locking and unlocking destroys performance. To address this limitation, most components pass batches of elements through the queues. As a component generates data, it saves it to an array of some fixed size. When the array is full, it is passed to the queue and a new array is created for the next batch. If the last batch is not full, it is copied to an array of an appropriate size so that there are not empty elements in the batch.

For a concrete example, consider the count module. It reads from source sequence files, k-merizes sequences, counts each instance of unique k-mers, and outputs the counts to a file. To count large data sets, a split component writes several "segments" of sorted k-mer counts to disk, and a merge component merges segments. This module loads and links the read, kmer, count split component, count merge component, and countfilewriter components. Data then flows from the source file to the destination file through this pipeline. The module itself does little more than create the components and link them

All components extend edu.gatech.kanalyze.comp.PipelineComponent. Component objects are created with their constructor and run with the parent class's run() method. Most components take two syncrhonized queues in their constructor (one for input, and one for output), and when run() is called, it begins reading from the input queue and writing to the output queue.

Input and output queues are parameterized instances of edu.gatech.kanalyze.util.BoundedQueue<T>. To signal the end of the stream, a null value is written to the stream. All components must re-write the null value back to the stream after reading it to ensure that all threads get the EOF signal.

The following sections outline the flow of information through each component. Any component specific properties, i.e. the one set with the -p command line option, are also documented.

#### 3.3.1 Count Merge Component

Input source: edu.gatech.kanalyze.util.BoundedQueue<File>

Output target: edu.gatech.kanalyze.util.BoundedQueue<CounterPair>

Class: edu.gatech.kanalyze.comp.count.CountMergeComponent

Complements the count split component (Section 3.3.2). Segment files created by the split component are merged and final k-mer counts are sent to the output queue.

Once null is read from the input queue, it is written back to the input queue (to ensure all threads get the signal).

Component properties:

**comp.count.deltemp** If set to the string "true" (not case sensitive), the segment files are deleted after they are read.

Global properties: None

## 3.3.2 Count Split Component

Input source: edu.gatech.kanalyze.util.BoundedQueue<long[]>
Output target: edu.gatech.kanalyze.util.BoundedQueue<File>
Class: edu.gatech.kanalyze.comp.count.CountSplitComponent

Complements the count merge component (Section 3.3.1). K-mers are received from the input queue and accumulated into a memory buffer until it is full. The memory buffer is sorted, and k-mer counts are written to a file on disk ("segment file"). The file location is sent on the output queue for the merge component. The memory buffer is cleared, and it receives the next batch of k-mers. The process is repeated until all k-mers are read.

Once null is read from the input queue, it is written back to the input queue (to ensure all threads get the signal).

Component properties:

comp.count.segsize Number of k-mers accumulated in the memory buffer before they are sorted and dumped to a segment file. This value must not be less than CountSplitComponent.MIN\_SEGMENT\_SIZE or greater than CountSplitComponent.MAX\_SEGMENT\_SIZE.

Global properties: None

#### 3.3.3 Count File Writer

Input source: edu.gatech.kanalyze.util.BoundedQueue<KmerCount>

Output target: String (File Name)

Class: edu.gatech.kanalyze.comp.countfilewriter.CountFileWriterComponent

Reads k-mer counts from an input queue and writes to a file with the specified file name.

Once null is read from the input queue, it is written back to the input queue (to ensure all threads get the signal), and the output file is closed.

Component properites: None

Global properties:

**kanalyze.outfmt** Controls the output format of k-mers. Values should be in, hex, or seq (default = seq).

#### 3.3.4 File Writer

Input source: edu.gatech.kanalyze.util.BoundedQueue<long[]>

Output target: String (File Name)

Class: edu.gatech.kanalyze.comp.filewriter.FileWriterComponent

Reads k-mers from the input queue and writes to a file. The file name is the output target. This component is useful for streaming k-mers into a file where counting or coalescing them is not required.

Once null is read from the input queue, it is written back to the input queue (to ensure all threads get the signal), and the output file is closed.

Component properites: None

Global properties:

**kanalyze.outfmt** Controls the output format of k-mers. Values should be in, hex, or seq (default = seq).

#### 3.3.5 K-mer

Input source: edu.gatech.kanalyze.util.BoundedQueue<SequenceRead[]>

Output target: edu.gatech.kanalyze.util.BoundedQueue<long[]>

Class: edu.gatech.kanalyze.comp.kmer.KmerComponent

Processes sequence reads from the input queue, k-merizes them, and writes k-mers to the output queue. Both the input queue and the output queue work in batches of items.

Once null is read from the input queue, it is written back to the input queue (to ensure all threads get the signal), and null is propagated to the output queue.

Component properites: None Global properties: None

#### 3.3.6 Reverse Complement

Input source: edu.gatech.kanalyze.util.BoundedQueue<long[]>
Output target: edu.gatech.kanalyze.util.BoundedQueue<long[]>

Class: edu.gatech.kanalyze.comp.rcompl.RevComplComponent

Reverse-complements k-mers from the input queue and writes to the output target. There are three modes (enum edu.gatech.comp.rcompl.RevComplMode), and one must be supplied to the constructor. If null, the default mode, RevComplMode.DUPLICATE is used.

Modes are as follows:

**RevComplMode.DUPLICATE** Outputs the original k-mer and its reverse complement. Most applications requiring reverse-complement will use this mode.

RevComplMode.LESSER Takes the reverse complement and outputs the lesser of the original or reverse complement (lesser defined by ASCII sort order of k-mers or the numeric sort order of numeric k-mers, which are the same). This mode may be useful for storing a set of k-mers and their reverse complements in about half the space required to store both.

**RevComplMode.SINGLE** Outputs the reverse complement of k-mers. The original k-mer is discarded.

Once null is read from the input queue, it is written back to the input queue (to ensure all threads get the signal), and null is propagated to the output queue.

Component properites: None Global properties: None

#### 3.3.7 Reader

Input source: edu.gatech.kmer.util.BoundedQueue<SequenceSource[]>

Output target: edu.gatech.kmer.util.BoundedQueue<SequenceRead[]>

Class: edu.gatech.kanalyze.comp.reader.ReaderComponent

Reads a queue of sequence sources, extracts sequence strings, and writes sequence read objects to the output queue.

Each SequenceSource object contains information about the sequence source (usually a file) and the declared type (e.g. FASTA or FASTQ). The type information is used to load the appropriate source reader. Source readers classes are located and objects instantiated using Java reflections. For more information about writing custom sequence readers for new file formats, see section 4.2.

Once null is read from the input queue, it is written back to the input queue (to ensure all threads get the signal), and null is propagated to the output queue.

Component properites: None

Global properties: None

# 3.4 Utility Classes

KAnalyze has several utility classes for shared functionality. Most of the classes are used by multiple components or modules. These classes may also be useful for other applications.

All of these classes are found in package edu.gatech.kanalyze.util

#### 3.4.1 BoundedQueue

BoundedQueue is a thread-safe queue of a fixed size. If the size is not set when the queue is created, the default size, edu.gatech.kanalyze.Constants.DEFAULT\_QUEUE\_SIZE is used. This class is parameterized, so it can queue objects of any type.

put() adds elements to the queue, and it will block until the queue is not full. take() gets elements from the queue, and it will block on an empty queue until an element is available.

Queues of this type are used to pass elements from one pipeline component to the next. It is designed to be simple and fast.

#### 3.4.2 KmerUtil

KmerUtil is a collection of k-mer methods used in several parts of the program. It performs tasks such as checking for a valid k size and converting k-mers among different representations.

#### 3.4.3 LongBoundedQueue

LongBoundedQueue was originally used for passing k-mers from one component to the next. Now that k-mers are passed in batches and batches are array objects, this class is no longer used in the main program. Since parts of the unit test code in edu.gatech.kanalyze.test still use it, and since it could be useful for some API, it has been left in the project. It may be removed in future revisions.

Besides passing long integers instead of classes, LongBoundedQueue is identical to BoundedQueue.

# 3.4.4 LongCounter

LongCounter is a hash-table based k-mer counter where the hash key is the k-mer and the value is the number of times the k-mer was added to the counter.

This class is very similar to java.util.HashTable in the standard Java API. In fact, the hash function is identical. However, because this implementation is not as general and because it places some restrictions on the bucket array (must be powers of 2), this implementation is generally twice as fast as Java HashTable.

This class is not currently used by any of the components packaged with KAnalyze, but it may be useful for future components or extensions outside of the main distribution.

# 3.4.5 LongHashSet

LongHashSet is a hash-table based set. It is used to keep track of which k-mers were encountered when the number of times each k-mer was found is not important.

Like LongCounter, this class is based on an optimized algorithm that outperforms implementation, java.util.HashSet.

This class is not currently used by any of the components packaged with KAnalyze, but it may be useful for future components or extensions outside of the main distribution.

## 3.4.6 LongIterator

LongIterator is an iterator interface that returns long integers instead of objects. For iterators over large k-mer sets, this iterator avoids the auto-boxing and auto-unboxing overhead of using the standard Java iterator interface. Normally, data structures that use this iterator also have an iterator method for retrieving a standard Java iterator.

This class is not currently used by any of the components packaged with KAnalyze, but it may be useful for future components or extensions outside of the main distribution.

## 3.4.7 StringUtil

StringUtil is a collection of string parsing utilities for parsing actions commonly found in KAnalyze. Currently, it only has the toNameCase() method.

toNameCase() is used by any class that uses a string to locate a class through the Java reflections API. For example, the sequence reader takes a string such as "fasta", which is not case sensitive, and uses it to find a class named "FastaReader". This method converts "fasta" to "Fasta".

# 4 Extending KAnalyze

KAnalyze was designed to be extensible. As part of this design, new features can be added to KAnalyze without modifying existing code. For example, to write a new module, sequence reader, or databse driver, no part of the existing system has know about the feature a priori.

First, consider the easy way to create the modular functionality of KAnalyze. From the command line, the first argument determines which module is loaded and executed. Typically, a programmer might write a main method that checks the first argument, uses it in a switch statement to match the string to the module name, and call the appropriate module. If a programmer wants to add a new module, the main method has to be updated so that it knows about the module.

KAnalyze takes a different approach. It does the work of finding modules on its own. This makes it easier to develop new modules, and it allows KAnalyze to find modules it doesn't know about at compile time. This means that any module can be run as long as the Java classloader can find it.

Modules are not the only part of KAnalyze that can do this. File readers can also be extended in this way to support new sequence file formats.

The key to this flexibility is the Java Reflection API. Given a name, such as "count", KAnalyze will find a module class, edu.gatech.kanalyze.module.count.CountModule. KAnalyze must be able to find the class and load it, so each of these extensible features has a set of requirements that must be met, such as what package the modules are in, what parent class they extend, and which constructors are available. The code running the feature must be able to instantiate and run the class polymorphically, hence the parent class and constructor requirements.

The Java class loader must be able to find classes loaded by the Java Reflection API. If two or more classes have the same name and are in the same package (perhaps in different JAR files), then it is up to the class loader to resolve the conflict.

#### 4.1 Modules

In the following set of requirements, xxx and Xxx is the name of the module in lower case and camel-case, respectively. Capitalization is imporant.

Module requirements:

- 1. Module class: edu.gatech.kanalyze.module.xxx.XxxModule
- 2. Module extends edu.gatech.kanalyze.module.KAnalyzeModule
- 3. Module has a constructor with the signature (java.lang.String[])
- 4. Module class is not abstract.

The KAnalyze default main method is edu.gatech.kanalyze.KAnalyzeModule.main(). This method reads its first argument and uses the class loader to find a class that adheres to the above rules. It is instantiated with the required constructor.

The main method then calls a series of methods. First, moduleInit(java.lang.String[]) accepts the command line arguments with the module name removed. The module's implementation of this method should use the arguments like any other main method. It should perform any required initialization tasks before returning. If moduleInit() does not throw an exception, moduleRun() is called.

After calling moduleInit() and moduleRun(), moduleFinalize() is called. If an exception was thrown by either of the previous two methods, it is passed to finalize so the finalize method in case it is important for the cleanup code to know if something failed. This finalize method is always called. It is also optional. The default finalize method does nothing, so module implementations may choose to override it. Note that unlike <code>java.lang.Class.finalize()</code>, this method is called as soon as the other two methods complete, not when the module object is garbage collected.

Module implementations should never call System.exit() directly or indirectly. If no errors or special conditions occur that stop the module from completing, the module should simply call return. If the module should terminate abnormally for any reason, it should throw an instance of edu.gatech.kanalyze.module.ModuleRunle.

Writing modules this way ensures that API users can call modules as part of larger programs. Therefore, custom module should *never* terminate the JVM with a call to System.exit().

If a help option is invoked on the command line, the module should print help information to standard output (java.lang.System.out) and throw a edu.gatech.kanalyze.module.ModuleRunException with return code edu.gatech.kanalyze.Constants.ERR\_NONE.

If any of the module's methods throws edu.gatech.kanalyze.module.ModuleRunException, main() will print an error to standard error, java.lang.System.err, message using the exception message and terminate the program with the return code in the exception. If the return code is edu.gatech.kanalyze.Constants.ERR\_main() assumes the program terminated normally (no error is output, and the program terminates with edu.gatech.kanalyze.Constants.ERR\_NONE).

# 4.2 Sequence Readers

In the following set of requirements, xxx and Xxx is the name of the driver in lower case and cammel-case, respectively. Capitalization is important.

Module requirements:

- 1. Reader class: edu.gatech.kanalyze.comp.reader.xxx.XxxSequenceReader
- 2. Reader extends edu.gatech.kanalyze.comp.reader.SequenceReader
- 3. Reader has a constructor with the signature (edu.gatech.kanalyze.util.BoundedQueue<edu.gatech.kanalyze.cint, Properties)
- 4. Reader class is not abstract.

In the following description, class names without a fully qualified package are found in package edu.gatech.kanalyze.com

The reader must implement the read() method from the parent class. Sequence Reader. The read method

The reader must implement the read() method from the parent class, SequenceReader. The read method takes a single argument, which is a SequenceSource object. The reader then calls getInputStream() from the source object and reads it until there is no more data to retrieve.

Each sequence string is then written to the next element of a batch array already setup by the parent class, readBatch[], while incrementing the batch counter, batchSize, until the batch counter is equal to the size of the batch array (the batch is full). The reader then calls flush(), which is implemented by the parent class. flush() will write the batch to the pipeline, create a new batch array, and reset the counter. To summarize, the reader writes sequences to a queue until it is full, then it calls flush, and it repeats the process until all sequences have been read.

The parent class has a protected variable, splitLength. If sequence reads are longer than this, they should be split into multiple reads. The last k - 1 elements from the end of the current split must be copied to the beginning of the next split to preserve correct k-mer generation. For fully assembled sequences, this prevents KAnalyze from memory starvation by attempting to input large reads at once.

Sequence readers should encapsulate the read() method in a try/finally that always calls flush() before it terminates. Failing to do so may result in loss of reads data.

# 5 Supplementary Information

# 5.1 Input File Formats

#### 5.1.1 FASTA

FASTA files are compatible with the FASTA format acceptable by NCBI BLAST<sup>4</sup>. The only exception is that the reader permits blank lines within the sequence because all blank lines are ignored by the parser.

A single-line description beginning with '¿' precedes each sequence record. The sequence record may be on a single line or split over any number of lines. There are no restrictions on the line length.

FASTA files are not checked by the reader. The sequence string may contain any characters as long a line does not begin with '¿'. The k-mer component that translates sequences skips over all unrecognized characters. Checking sequences while reading is omitted for performance reasons.

Standard line break characters are recognized: LF (Linux/Unix/MacOS), or CRLF (Windows). The file may or may not end with a newline character.

#### 5.1.2 FASTAGZ

A GZIP compressed FASTA file. The uncompressed file contents must follow the FASTA format defined in Section 5.1.1.

#### 5.1.3 FASTQ

FASTQ files must be compatible with the format published in Nucleic Acid Reseach<sup>5</sup>. The first line of a record begins with '@'. The sequence may be split over multiple lines. A line beginning with '+' will separate the sequence from the quality scores. Quality scores are read line-by-line until the number of quality score characters is the same as the number of sequence characters read for any given record. The quality score characters are discarded (this reader does not enforce proper scores).

This reader permits blank lines between records, but not within records.

Standard line break characters are recognized: LF (Linux/Unix/MacOS), or CRLF (Windows). The file may or may not end with a newline character.

#### 5.1.4 FASTQGZ

A GZIP compressed FASTQ file. The uncompressed file contents must follow the FASTQ format defined in Section 5.1.3.

## 5.1.5 RAW

Raw sequence files are the simplest record type. A sequence may be split over any number of lines. Blank lines separate records.

Standard line break characters are recognized: LF (Linux/Unix/MacOS), or CRLF (Windows). The file may or may not end with a newline character.

<sup>&</sup>lt;sup>4</sup>http://www.ncbi.nlm.nih.gov/BLAST/blastcgihelp.shtml

<sup>&</sup>lt;sup>5</sup>http://www.ncbi.nlm.nih.gov/pmc/articles/PMC2847217/

# 5.2 Kmer Output Formats

K-mers may be output in one of several formats. The most natural is the sequence format, which a sequence of k bases, each A, C, T, or G. Note that if RNA sequences are read, U is converted to T on output. This is because T and U are represented the same internally.

Two other representations, decimal and hexadecimal, take a little more explanation. Internally, k-mers are represented as integers. Converting the sequences to integers makes processing and storing the k-mers fast and efficient.

To understand the k-mer nucletide sequence to integer conversion, the nucleotide sequence can be thought of a base-4 numbering system where A is 0, C is 1, G is 2, and T is 3. Each base can be conveniently represented as a two-bit number where A = 00, C = 01, G = 10, and G = 11. A k-mer sequence is then a sequence of these letters. For example, the 8-mer CCAGTCGT is 0101001011011011.

Note that k-mers are given a natural order where lower numbers come before higher numbers. For 8-mers, the first possible k-mer is AAAAAAAA, and the last possible k-mer is TTTTTTT. This numeric ordering can be taken advantage of in many ways.

Java's long integer is a 64-bit 2's complement integer. Each k-mer is is two bits in length. In order to preserve the natural ordering, we chose not to include negative numbers, which is any 2's complement integer that begins with 1. Excluding the first bit, which is always 0, a 64-bit integer can represent any k-mer size 1 to 31. Therefore, the maximum k-mer size is currently 31.

Binary numbers are conveniently represented in hexadecimal. K-mers can also be written as hexadecimal integers, which affords them a compact representation. For example, the same 8-mer as above, CCAGTCGT, can be represented as 0x52DB.

The global property, *kanalyze.outfmt*, sets one of the following output options. If the property is not set, the default sequence option is set. For convenience, most modules have a shorter command line option to set this property.

# 5.2.1 Sequence Format

Tab delimited file of k-mers and counts. K-mers are written as a string of A, C, G, and T. Each k-mer has k letters. Note that if RNA sequences were read, U will be converted to T when sequences are output in this format.

#### 5.2.2 Decimal Format

Tab delimited file of k-mers and counts using the base-10 decimal representation of k-mers. For more information about how k-mer sequences are represented as integers, see the above section 5.2.

## 5.2.3 Hexadecimal Format

Tab delimited file of k-mers and counts using the base-16 decimal representation of k-mers. For more information about how k-mer sequences are represented as integers, see the above section 5.2.

## 5.2.4 FASTA Format

FASTA representation of k-mers counts. Each FASTA record is a single k-mer and its count. The description line is the k-mer count and the sequence is the k-mer sequence.

#### 5.3 Command Line Return Codes

The KAnalyze user interface always returns a well-defined code. When executed from a script environment, this return code can easily be checked to see if KAnalyze completed normally or not. Each return code is defined in this section.

In the following list of return codes, the numeric code is listed. For API users, the constant defined in edu.gatech.kanalyze.Constants is also listed.

- **0** (ERR\_NONE) The program terminated normally. All k-mers were successfully processed without error, or the help option was invoked.
- 1 (ERR\_USAGE) Command line arguments were incomplete, improperly formatted, or required arguments were missing.
- **2 (ERR\_IO)** An I/O (input/output) error occured reading or writing data. This error is normally returned for file I/O errors.
- **3 ERR\_SECURITY** A security error, such as permissions denied, occured.
- 4 (ERR\_FILENOTFOUND) A required or specified file was not found.
- **5** (ERR\_DATAFORMAT) Data in an input file is improperly formatted.
- 6 (ERR\_ANALYSIS) Some un-recoverable error occured during analysis.
- 7 (ERR\_INTERRUPTED) A program thread was interrupted while it was running. This would normally be returned if the program is terminated before it completes.
- 99 (ERR\_SYSTEM) Some serious unrecoverable system error occured. These errors are almost certainly program bugs. Some other unrecoverable errors, such as running out of memory, could return this code.

# 5.4 Building From Source

The source is distributed with an Apache Ant build file with multiple targets. To build these, Apache Ant must be installed.

To run a target, enter ant <target>. Multiple targets may be supplied, for example ant clean package.

The most commonly used targets are "clean", "compile", and "package". The clean removes all temporary files. It is not often necessary, but it should be executed before building a package for testing to ensure no artifacts are left behind from previous builds. The compile target compiles the class files and does nothing else. The package target generates the JAR file and distributable packages.

The "doc.javadoc" target generates the Javadoc pages from KAnalyze source comments. See Section 5.4.1 for more information about these pages. "doc.manual" generates this manual as a PDF file.

The "test" target executes unit tests packaged with KAnalyze. See Section 5.5 for more information on unit testing.

Other targets are called by the targets described above and are not often called directly.

For a full list of targets and brief descriptions, enter ant -projecthelp.

#### 5.4.1 Building Javadoc Pages

Javadoc is a very powerful tool for documenting APIs written in Java. Section 3.1 introduces the concept and the KAnalyze rules governing these comments.

Javadoc comments begin with "/\*\*", and end with "\*/". In KAnalyze, these comments appear before every class, method, and field. The comment starts with a description of what the element does. For methods, it documents parameters and the conditions under which all exceptions are thrown. Full documentation for writing these comments can be found online<sup>6</sup>.

The ant build system (Section 5.4) generates web pages from these comments. The "doc.javadoc" target builds two sets of pages. One, the API documentation, documents all elements available on the API. This target is intended for developers who are extending or using the API. The other, FULL documentation, documents everything including private members that are not available to the API. This is meant for KAnalyze maintenance programmers.

 $<sup>^6 \</sup>rm http://www.oracle.com/technetwork/java/javase/documentation/index-137868.html$ 

After running the javadoc Ant target, the API documentation can be access by loading "build/doc/javadoc/api/index.html in a browser (relative to the project root). The FULL documentation can be access by loading "build/doc/javadoc/full/index in a browser (relative to the project root).

# 5.5 Running Unit Tests

KAnalyze is distributed with a set of unit tests. As changes are made to the system, existing components can be quickly re-tested. These tests are intended to be thorough. For example, the unit tests for the k-mer component tests all valid k-mer sizes. Some of the utility class tests use up to 8GB of memory.

Tests are easily executed with "ant test". This target compiles the test classes and executes all of them. Note that test classes are not compiled by the "compile" target and are not distributed with compiled packages.

# 6 License

# 6.1 Documetation

This document is licensed under the GNU Free Documentation License 1.3 or later. The file "COPY-ING.DOC" contains the text of this license. If you did not receive the file in the source distribution, please contact us or the GNU website for a copy.

# 6.2 KAnalyze Software

The KAnalyze software is licensed under the GNU Lesser General Public License version 3 or later. The file "COPYING" contains the text of the GNU GPL, and the file "COPYING.LESSER" contains the GNU LGPL extension to the GPL. If you did not receive the file in the source distribution, please contact us or the GNU website for a copy.